GOLDEN SEASON PTE. LTD.

www.goldenseason.com.sg | www.gseason.com



Body Bags



The Body Bag provides a means to store corpses and remains, as well as transportation for post-mortem or clearing of disaster sites. This is a one size fits all body bag.

Water proof material

Golden Season's Body bags are manufactured using high quality PVC ensuring the water proofing and no leaking of fluids, which may cause contamination.

Flexible, Compact and Easy to store and transport

The body bags are foldable and easily deployed for usage.

Heavy-Duty Zippers

The body bags are equipped with heavy-duty zippers to allow access to the body for inspection. The body does not have to be removed from the bag for inspection and post-mortem.

Carrying Handles

The body bag comes with heavy-duty handles to aid transportation.

| Specifications | |
|-------------------|---------------------------------------|
| Material: | Heavy Duty PVC |
| Carrying Handles: | Yes, 6 Handles |
| Dimension: | 2200mm x 900mm x 5mm (Empty) |
| Accessories: | Carrying handles (6 provided per set) |

^{*}Specifications are subject to change without prior notice.