GOLDEN SEASON PTE. LTD.



Civilian Stretcher Stand







The Civilian stretcher stand is used in conjunction with our other stretcher to elevate the resting height of patients and casualties. It assists paramedics and health care professionals to examine the patients at waist level versus low ground level.

Compactible

The litter stand is compactible for easy storage.

Easy to operate

Golden Season's civilian stretcher stand is easy to operate.

Lightweight

Golden Season's Civilian Stretcher Stand only 2.21 kg. Each.

| Specifications | |
|-------------------|---|
| Dimension Open: | 61cm x 48cm x 87cm |
| Dimension Folded: | 61cm x 7cm x 90cm |
| Weight: | 2.21kg each litter stand (approx. 4.43kg each pair) |
| Material Type: | Lightweight Aluminium |
| Accessories: | 2 folding Stretcher (Sold Separately) |

^{*}Specifications are subject to change without prior notice.