GOLDEN SEASON PTE. LTD.



Flexible Splint









The flexible splint is moldable and allows for shaping into a strong support for any fractured or injured limb, as well as neck injuries.

Compact

The Flexible splint rolls up into a compact bundle for easy storage and transportation.

X-ray Translucent

Manufactured from food grade alloy, the flexible splint allows for X-Ray without removal.

Easily shaped and cut to size

Golden Season's Flexible Splint can be cut and shaped to any form to fit the different fracture locations.

Specifications	
Width:	90cm
Height:	11cm
Weight:	145g
Accessories:	Bandages(Sold Separately)

^{*}Specifications are subject to change without prior notice.