



ACE Cervical Collar



The Perfit ACE® Cervical Collar is constructed of precision-engineered, heavy-duty plastic and lined with a protective layer of foam to ensure both immobilization and comfort. The back of the collar is elongated to give complete support. However, the feature that makes it the neck support of choice is its adaptability. With 16 different settings, this cervical collar provides a perfect fit for any adult casualty.

Special Features:

- ✓ 16 adjustable settings
- ✓ Foamed collar
- ✓ Tactical or white color

Contents:

Item#: 50-0010 NSN#: 6515-01-541-8147 (Tactical Color)

Item#: 50-1007 (White Color)

1 x ACE Cervical Collar

Specifications	
Dimensions:	H 6.75 in. x W 20.5 in. x D 0.75 in
Weight:	5.2 oz

*Specifications are subject to change without prior notice.

Golden Season Pte Ltd