



NAR Hypothermia Prevention and Management Kit (HPMK)



NOTE: Litter Not Included

Recommended by CoTCCC for the prevention of hypothermia during casualty care, the patented HPMK® with reinforced Heat Reflective Shell is strong, flexible, lightweight & impervious to wind & rain. The HRS™, the kit's ruggedized outer shell, is constructed of a 4-ply, composite fabric with a protected non-conductive reflective layer that provides excellent thermal insulation. Rapid 360° access to the casualty with minimal exposure to the elements is achieved through the use of 1.5 in. continuous Velcro® closures located along the perimeter of the covering.

The HRS™ also features a built-in hood, fluid absorption pad & top-to-bottom tapered shape to maximize isothermal capabilities. The heat-generating source of the HPMK® is a self-heating, oxygen-activated shell liner designed to sustain 6 hours of continuous dry heat at 106°F (no external power supply required). This low cube kit is vacuum-packaged in a rugged, pliable pouch featuring NAR's Red Tip Technology® with tabs that allow easy access and rapid deployment.

Product of North American Rescue



Golden Season Pte Ltd



Special Features:

- ✓ TCCC-recommended hypothermia prevention solution during casualty care
- ✓ Sustains 6 hours of continuous dry heat to 106°F with oxygen-activated, self-heating liner (no external power supply required)
- ✓ Constructed of 4-ply composite fabric with protected non-conductive reflective layer that provides excellent thermal insulation
- ✓ Tapered shape from top to bottom to maximize isothermal capabilities
- ✓ Built-in hood & integrated fluid absorption pad
- ✓ Continuous 1.5 in. Velcro® closures for rapid 360° access to the casualty with minimal exposure to the elements
- ✓ Reinforced and weight tested to support 250+ lb
- ✓ Vacuum-packaged in a durable pouch for a protective, low-cube footprint
- ✓ Storage pouch features NAR's Red Tip Technology® with tabs that allow easy access and rapid deployment
- ✓ U.S. Patent 7,766,950 B2

Contents:

- ✓ 1 x Heat Reflective Shell (HRS™) with built-in hood
- ✓ 1 x Self-heating shell liner
- ✓ 1 x Plastic vacuum bag

Specifications	
Packaged (HPMK®):	H 6.75 in. x W 10.5 in. x D 5.5 in.
Weight (HPMK®):	3 lb 8 oz
Heat Reflective Shell (Open):	L 78 in. x W 43 in.

*Specifications are subject to change without prior notice.